

For Research Use Only

# NLRP3 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68102-3-PBS

## Basic Information

<b>Catalog Number:</b> 68102-3-PBS	<b>GenBank Accession Number:</b> NM_001079821	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 114548	<b>CloneNo.:</b> 2F1B11
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q96P20	
<b>Isotype:</b> IgG1	<b>Full Name:</b> NLR family, pyrin domain containing 3	
<b>Immunogen Catalog Number:</b> AG26289	<b>Calculated MW:</b> 118 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP80005-1, NLRP3 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 80422-2-PBS. Detection antibody: 68102-3-PBS. Standard: Ag26289. Range: 1.563-25 ng/mL