

For Research Use Only

# SLC25A16 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68107-4-PBS



## Basic Information

Catalog Number:

68107-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4891

GenBank Accession Number:

BC030266

GeneID (NCBI):

8034

UNIPROT ID:

P16260

Full Name:

solute carrier family 25  
(mitochondrial carrier; Graves disease  
autoantigen), member 16

Calculated MW:

36 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G9C7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50631-2, SLC25A16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68107-4-PBS. Detection antibody: 68107-5-PBS. Standard:Ag4891. Range: 0.195-100 ng/mL



Cytometric bead array standard curve of MP50631-3, SLC25A16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68107-4-PBS. Detection antibody: 68107-6-PBS. Standard:Ag4891. Range: 0.781-100 ng/mL