

For Research Use Only

SLC25A16 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68107-6-PBS



Basic Information

Catalog Number: 68107-6-PBS	GenBank Accession Number: BC030266	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 8034	CloneNo.: 3F9B5
Source: Mouse	UNIPROT ID: P16260	
Isotype: IgM	Full Name: solute carrier family 25 (mitochondrial carrier; Graves disease autoantigen), member 16	
Immunogen Catalog Number: AG4891	Calculated MW: 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50631-3, SLC25A16 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68107-4-PBS. Detection antibody: 68107-6-PBS. Standard: Ag4891. Range: 0.781-100 ng/mL