

For Research Use Only

AXL Monoclonal antibody, PBS Only

Catalog Number: 68108-1-PBS



Basic Information

Catalog Number: 68108-1-PBS	GenBank Accession Number: BC032229	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 558	CloneNo.: 3D12F11
Source: Mouse	UNIPROT ID: P30530	
Isotype: IgG1	Full Name: AXL receptor tyrosine kinase	
Immunogen Catalog Number: AG14401	Calculated MW: 894 aa, 98 kDa	
	Observed MW: 130-140 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68108-1-Ig (AXL antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with Lamin B1 Monoclonal antibody (66095-1-Ig) as loading control. This data was developed using the same antibody clone with 68108-1-PBS in a different storage buffer formulation.