

For Research Use Only

PGAM5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68116-3-PBS



Basic Information

Catalog Number:

68116-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28195

GenBank Accession Number:

BC008196

GeneID (NCBI):

192111

UNIPROT ID:

Q96HS1

Full Name:

phosphoglycerate mutase family
member 5

Calculated MW:

32 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B7B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50580-1, PGAM5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68116-2-PBS. Detection antibody: 68116-3-PBS. Standard:Ag28195. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50580-2, PGAM5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68116-4-PBS. Detection antibody: 68116-3-PBS. Standard:Ag28195. Range: 1.563-100 ng/mL.