

For Research Use Only

PGAM5 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68116-6-PBS



Basic Information

| | | |
|---|---|--|
| Catalog Number: 68116-6-PBS | GenBank Accession Number: BC008196 | Purification Method: Protein G Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 192111 | CloneNo.: 2E10F1 |
| Source: Mouse | UNIPROT ID: Q96HS1 | |
| Isotype: IgG1 | Full Name: phosphoglycerate mutase family member 5 | |
| Immunogen Catalog Number: AG28195 | Calculated MW: 32 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50580-3, PGAM5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68116-5-PBS. Detection antibody: 68116-6-PBS. Standard:Ag28195. Range: 1.563-100 ng/mL.