

For Research Use Only

Asfv P30 of ASFV Monoclonal antibody, PBS Only



Catalog Number: 68134-1-PBS

Basic Information

Catalog Number: 68134-1-PBS	GenBank Accession Number: M96354	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 22220322	CloneNo.: 3B9E3
Source: Mouse	UNIPROT ID: P34204	
Isotype: IgG1	Full Name: Phosphoprotein p30	
Immunogen Catalog Number: AG31446	Calculated MW: 30 kDa	
	Observed MW: 35 kDa	

Applications

Tested Applications:
WB, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein of P30 of ASFV This data was developed using the same antibody clone with 68134-1-PBS in a different storage buffer formulation.