

For Research Use Only

NDUFV2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68150-2-PBS



Basic Information

Catalog Number:

68150-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7760

GenBank Accession Number:

BC001632

GeneID (NCBI):

4729

UNIPROT ID:

P19404

Full Name:

NADH dehydrogenase (ubiquinone)
flavoprotein 2, 24kDa

Calculated MW:

27 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F9A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50478-1, NDUFV2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68150-2-PBS. Detection antibody: 68150-3-PBS. Standard:Ag7760. Range: 0.195-25 ng/mL



Cytometric bead array standard curve of MP50478-2, NDUFV2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68150-2-PBS. Detection antibody: 68150-4-PBS. Standard:Ag7760. Range: 0.391-00 ng/mL