

For Research Use Only

# DDAH2 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68154-3-PBS



## Basic Information

Catalog Number: 68154-3-PBS	GenBank Accession Number: BC001435	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 23564	CloneNo.: 2C3A8
Source: Mouse	UNIPROT ID: O95865	
Isotype: IgG1	Full Name: dimethylarginine dimethylaminohydrolase 2	
Immunogen Catalog Number: AG7127	Calculated MW: 29 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50417-1, DDAH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68154-2-PBS. Detection antibody: 68154-3-PBS. Standard:Ag7127. Range: 0.195-100 ng/mL.



Cytometric bead array standard curve of MP50417-3, DDAH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68154-3-PBS. Detection antibody: 68154-4-PBS. Standard:Ag7127. Range: 0.195-100 ng/mL.