

For Research Use Only

DDAH2 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 68154-4-PBS



Basic Information

Catalog Number: 68154-4-PBS	GenBank Accession Number: BC001435	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 23564	CloneNo.: 1G8G9
Source: Mouse	UNIPROT ID: O95865	
Isotype: IgG1	Full Name: dimethylarginine dimethylaminohydrolase 2	
Immunogen Catalog Number: AG7127	Calculated MW: 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50417-2, DDAH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68154-4-PBS. Detection antibody: 68154-1-PBS. Standard:Ag7127. Range: 0.391-50 ng/mL



Cytometric bead array standard curve of MP50417-3, DDAH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68154-3-PBS. Detection antibody: 68154-4-PBS. Standard:Ag7127. Range: 0.195-100 ng/mL