

For Research Use Only

HIBCH Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68157-5-PBS



Basic Information

Catalog Number:

68157-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6473

GenBank Accession Number:

BC067822

GeneID (NCBI):

26275

UNIPROT ID:

Q6NVY1

Full Name:

3-hydroxyisobutyryl-Coenzyme A
hydrolase

Calculated MW:

43 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C7E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50627-3, HBCH Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68157-5-PBS. Detection antibody: 68157-3-PBS. Standard: Ag6473. Range: 0.195-100 ng/mL.