

For Research Use Only

ERBB2IP Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68159-2-PBS

Basic Information

Catalog Number: 68159-2-PBS	GenBank Accession Number: BC115012	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55914	CloneNo.: 1C12G9
Source: Mouse	UNIPROT ID: Q96RT1	
Isotype: IgG1	Full Name: erbb2 interacting protein	
Immunogen Catalog Number: AG18514	Calculated MW: 1412 aa, 158 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51110-1, ERBB2IP Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68159-1-PBS. Detection antibody: 68159-2-PBS. Standard:Ag18514. Range: 0.195-100 ng/mL.