

For Research Use Only

GPD2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68174-2-PBS



Basic Information

Catalog Number: 68174-2-PBS	GenBank Accession Number: BC019874	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 2820	CloneNo.: 1E10D6
Source: Mouse	UNIPROT ID: P43304	
Isotype: IgG2a	Full Name: glycerol-3-phosphate dehydrogenase 2 (mitochondrial)	
Immunogen Catalog Number: AG11212	Calculated MW: 41 kDa, 81 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50432-1, GPD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68174-2-PBS. Detection antibody: 68174-1-PBS. Standard:Ag11212. Range: 0.391-200 ng/mL



Cytometric bead array standard curve of MP50432-2, GPD2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68174-2-PBS. Detection antibody: 68174-3-PBS. Standard:Ag11212. Range: 0.781-200 ng/mL