

For Research Use Only

DFFA,DFF45 Monoclonal antibody, PBS Only (Capture)



Catalog Number:68177-3-PBS

Basic Information

| | | |
|---|--|---|
| Catalog Number: 68177-3-PBS | GenBank Accession Number: BC000037 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 1676 | CloneNo.: 2B10B11 |
| Source: Mouse | UNIPROT ID: O00273 | |
| Isotype: IgG1 | Full Name: DNA fragmentation factor, 45kDa, alpha polypeptide | |
| Immunogen Catalog Number: AG32657 | Calculated MW: 37 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50354-2, DFFA, DFF45 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68177-3-PBS. Detection antibody: 68177-1-PBS. Standard: Ag32657. Range: 6.25-200 ng/mL.