

For Research Use Only

# PSMD7 Monoclonal antibody, PBS Only



Catalog Number: 68188-1-PBS

Featured Product

## Basic Information

Catalog Number:

68188-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9305

GenBank Accession Number:

BC012606

GeneID (NCBI):

5713

UNIPROT ID:

P51665

Full Name:

proteasome (prosome, macropain)

26S subunit, non-ATPase, 7

Calculated MW:

324 aa, 37 kDa

Observed MW:

37 kDa

Purification Method:

Protein G purification

CloneNo.:

3E10E7

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human

## Background Information

Proteasome 26S subunit, non-ATPase 7 (PSMD7), an ATP-independent component of the 19S regulatory subunit, is a member of the JAMM/MPN domain-associated metallopeptidase (JAMM) DUB family. PSMD7, as a core component of the 26S proteasome, is critical for the degradation of ubiquitinated proteins in the proteasome. PSMD7 arrests cell cycle in the G2/M phase during HIV infection. The molecular weight of PSMD7 is 37 kDa. (PMID: 34512150, 34234864)

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68188-1-Ig (PSMD7 antibody) at dilution of 1:20000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 68188-1-PBS in a different storage buffer formulation.



WB result of PSMD7 antibody (68188-1-Ig; 1:30000; incubated at room temperature for 1.5 hours) with sh-Control and sh-PSMD7 transfected HeLa cells. This data was developed using the same antibody clone with 68188-1-PBS in a different storage buffer formulation.