

For Research Use Only

# NNT Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68194-4-PBS



## Basic Information

Catalog Number:

68194-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4259

GenBank Accession Number:

BC032370

GeneID (NCBI):

23530

UNIPROT ID:

Q13423

Full Name:

nicotinamide nucleotide  
transhydrogenase

Calculated MW:

1085 aa, 114 kDa

Purification Method:

Protein G purification

CloneNo.:

3H5E5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50592-2, NNT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68194-4-PBS. Detection antibody: 68194-3-PBS. Standard:Ag4259. Range: 0.195-12.5 ng/mL