

For Research Use Only

IDH3B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68199-4-PBS



Basic Information

Catalog Number: 68199-4-PBS	GenBank Accession Number: BC001960	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 3420	CloneNo.: 1H3F12
Source: Mouse	UNIPROT ID: O43837	
Isotype: IgG1	Full Name: isocitrate dehydrogenase 3 (NAD+) beta	
Immunogen Catalog Number: AG8693	Calculated MW: 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50690-2, IDH3B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68199-2-PBS. Detection antibody: 68199-4-PBS. Standard: Ag8693. Range: 3.125-200 ng/mL.