

For Research Use Only

# Integrin alpha-1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68211-4-PBS



## Basic Information

Catalog Number:

68211-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30631

GenBank Accession Number:

BC137121

GeneID (NCBI):

3672

UNIPROT ID:

P56199

Full Name:

integrin, alpha 1

Calculated MW:

1179 aa, 131 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2F12E9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50423-2, Integrin alpha-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68211-4-PBS. Detection antibody: 68211-3-PBS. Standard: Ag30631. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50423-3, Integrin alpha-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68211-4-PBS. Detection antibody: 68211-1-PBS. Standard: Ag30631. Range: 0.781-100 ng/mL.