

For Research Use Only

# NF90,ILF3 Monoclonal antibody, PBS Only (Detector)

Catalog Number:68213-4-PBS



## Basic Information

Catalog Number:

68213-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30868

GenBank Accession Number:

BC064836

GeneID (NCBI):

3609

UNIPROT ID:

Q12906

Full Name:

interleukin enhancer binding factor 3,

90kDa

Calculated MW:

894 aa, 95 kDa

Purification Method:

Protein G purification

CloneNo.:

1A9E3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50472-3, NF90,ILF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68213-1-PBS. Detection antibody: 68213-4-PBS. Standard:Ag30868. Range: 0.781-100 ng/mL