

For Research Use Only

NF90,ILF3 Monoclonal antibody, PBS Only (Detector)

Catalog Number:68213-6-PBS



Basic Information

Catalog Number: 68213-6-PBS	GenBank Accession Number: BC064836	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 3609	CloneNo.: 2G6C2
Source: Mouse	UNIPROT ID: Q12906	
Isotype: IgG1	Full Name: interleukin enhancer binding factor 3, 90kDa	
Immunogen Catalog Number: AG30868	Calculated MW: 894 aa, 95 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50472-4, NF90,ILF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68213-5-PBS. Detection antibody: 68213-6-PBS. Standard:Ag30868. Range: 1563-100 ng/mL