

For Research Use Only

NF90,ILF3 Monoclonal antibody, PBS Only (Capture)

Catalog Number:68213-7-PBS



Basic Information

Catalog Number:

68213-7-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30868

GenBank Accession Number:

BC064836

GeneID (NCBI):

3609

UNIPROT ID:

Q12906

Full Name:

interleukin enhancer binding factor 3,

90kDa

Calculated MW:

894 aa, 95 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A3A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50472-2, NF90,ILF3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68213-7-PBS. Detection antibody: 68213-3-PBS. Standard:Ag30868. Range: 0.781-50 ng/mL