

For Research Use Only

RTN3 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68215-2-PBS

Basic Information

Catalog Number:

68215-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG13766

GenBank Accession Number:

BC011394

GeneID (NCBI):

10313

UNIPROT ID:

O95197

Full Name:

reticulon 3

Calculated MW:

113 kDa

Purification Method:

Protein A purification

CloneNo.:

1D7F9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50025-1, RTN3
Monoclonal Matched Antibody Pair, PBS Only.
Capture antibody: 68215-2-PBS. Detection
antibody: HRP-conjugated 68215-1-PBS. Standard:
Ag13766. Range: 0.625-20 ng/mL.