

For Research Use Only

HIGD1A Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68231-4-PBS



Basic Information

Catalog Number:

68231-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14027

GenBank Accession Number:

BC000601

GeneID (NCBI):

25994

UNIPROT ID:

Q9Y241

Full Name:

HIG1 hypoxia inducible domain
family, member 1A

Calculated MW:

93 aa, 10 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1A7A6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50454-2, HIGD1A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68231-1-PBS. Detection antibody: 68231-4-PBS. Standard: Ag14027. Range: 0.391-100 ng/mL.