

For Research Use Only

GPSM1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68233-2-PBS



Basic Information

Catalog Number:

68233-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG18773

GenBank Accession Number:

BC017353

GeneID (NCBI):

26086

UNIPROT ID:

Q86YR5

Full Name:

G-protein signaling modulator 1
(AGS3-like, C. elegans)

Calculated MW:

652 aa, 72 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C9C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50909-1, GPSM1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68233-2-PBS. Detection antibody: 68233-3-PBS. Standard: Ag18773. Range: 12.5-200 ng/mL