

For Research Use Only

# ALDH2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68237-3-PBS



## Basic Information

Catalog Number:

68237-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG3

Immunogen Catalog Number:

AG7463

GenBank Accession Number:

BC002967

GeneID (NCBI):

217

UNIPROT ID:

P05091

Full Name:

aldehyde dehydrogenase 2 family  
(mitochondrial)

Calculated MW:

56 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B5A3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50663-1, ALDH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68237-2-PBS. Detection antibody: 68237-3-PBS. Standard: Ag7463. Range: 0.781-200 ng/mL.