

For Research Use Only

ALDH2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68237-4-PBS



Basic Information

Catalog Number:

68237-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG7463

GenBank Accession Number:

BC002967

GeneID (NCBI):

217

UNIPROT ID:

P05091

Full Name:

aldehyde dehydrogenase 2 family
(mitochondrial)

Calculated MW:

56 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2E12G4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50663-2, ALDH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68237-4-PBS. Detection antibody: 68237-5-PBS. Standard: Ag7463. Range: 3.125-200 ng/mL.



Cytometric bead array standard curve of MP50663-4, ALDH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68237-4-PBS. Detection antibody: 68237-8-PBS. Standard: Ag7463. Range: 1.563-200 ng/mL.