

For Research Use Only

ALDH2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68237-7-PBS



Basic Information

Catalog Number: 68237-7-PBS	GenBank Accession Number: BC002967	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 217	CloneNo.: 2H8B10
Source: Mouse	UNIPROT ID: P05091	
Isotype: IgG1	Full Name: aldehyde dehydrogenase 2 family (mitochondrial)	
Immunogen Catalog Number: AG7463	Calculated MW: 56 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50663-3, ALDH2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68237-6-PBS. Detection antibody: 68237-7-PBS. Standard: Ag7463. Range: 3.125-200 ng/mL.