

For Research Use Only

PDHB Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68238-3-PBS



Basic Information

Catalog Number:

68238-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6857

GenBank Accession Number:

BC000439

GeneID (NCBI):

5162

UNIPROT ID:

P11177

Full Name:

pyruvate dehydrogenase (lipoamide)

beta

Calculated MW:

39 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1E4H8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50446-1, PDHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68238-2-PBS. Detection antibody: 68238-3-PBS. Standard:Ag6857. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50446-2, PDHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68238-4-PBS. Detection antibody: 68238-3-PBS. Standard:Ag6857. Range: 0.391-100 ng/mL.