

For Research Use Only

FBXO17 Monoclonal antibody, PBS Only



Catalog Number: 68246-1-PBS

Basic Information

Catalog Number: 68246-1-PBS	GenBank Accession Number: BC012385	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 115290	CloneNo.: 1F3D3
Source: Mouse	UNIPROT ID: Q96EF6	
Isotype: IgG1	Full Name: F-box protein 17	
Immunogen Catalog Number: AG3867	Calculated MW: 278 aa, 31 kDa	
	Observed MW: 29-31 kDa	

Applications

Tested Applications:
WB, Indirect ELISA, IF

Species Specificity:
Human

Background Information

FBXO17, also named as FBG4, FBX17, FBX26 and FBXO26, is a substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



NCI-H1299 cells were subjected to SDS PAGE followed by western blot with 68246-1-Ig (FBXO17 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68246-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using FBXO17 antibody (68246-1-Ig, Clone: 1F3D3) at dilution of 1:800 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L), CL594-Phalloidin (red). This data was developed using the same antibody clone with 68246-1-PBS in a different storage buffer formulation.