

For Research Use Only

TRA-1-60 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68250-2-PBS



Basic Information

Catalog Number:

68250-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12859

GenBank Accession Number:

BC093730

GeneID (NCBI):

5420

UNIPROT ID:

O00592

Full Name:

podocalyxin-like

Calculated MW:

526 aa, 55 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2A2A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51166-1, TRA-1-60 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68250-2-PBS. Detection antibody: 68250-3-PBS. Standard:Ag12859. Range: 0.391-25 ng/mL



Cytometric bead array standard curve of MP51166-2, TRA-1-60 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68250-2-PBS. Detection antibody: 68250-4-PBS. Standard:Ag12859. Range: 0.391-100 ng/mL