

For Research Use Only

TRA-1-60 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68250-4-PBS



Basic Information

Catalog Number: 68250-4-PBS	GenBank Accession Number: BC093730	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 5420	CloneNo.: 4H10H4
Source: Mouse	UNIPROT ID: O00592	
Isotype: IgG1	Full Name: podocalyxin-like	
Immunogen Catalog Number: AG12859	Calculated MW: 526 aa, 55 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51166-2, TRA-1-60 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68250-2-PBS. Detection antibody: 68250-4-PBS. Standard:Ag12859. Range: 0.391-100 ng/mL.