

For Research Use Only

NDUFS1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68253-5-PBS



Basic Information

Catalog Number:

68253-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG3135

GenBank Accession Number:

BC030833

GeneID (NCBI):

4719

UNIPROT ID:

P28331

Full Name:

NADH dehydrogenase (ubiquinone)

Fe-S protein 1, 75kDa (NADH-coenzyme Q reductase)

Calculated MW:

727 aa, 79 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3C11G6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50476-2, NDUFS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68253-4-PBS. Detection antibody: 68253-5-PBS. Standard: Ag3135. Range: 0.098-12.5 ng/mL