

For Research Use Only

NDUFS1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68253-6-PBS

Basic Information

Catalog Number: 68253-6-PBS	GenBank Accession Number: BC030833	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4719	CloneNo.: 2A7B1
Source: Mouse	UNIPROT ID: P28331	
Isotype: IgG2b	Full Name: NADH dehydrogenase (ubiquinone) Fe-S protein 1, 75kDa (NADH- coenzyme Q reductase)	
Immunogen Catalog Number: AG3135	Calculated MW: 727 aa, 79 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50476-3, NDUFS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68253-6-PBS. Detection antibody: 68253-7-PBS. Standard: Ag3135. Range: 0.098-3.125 ng/mL.