

For Research Use Only

# NDUFS1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68253-8-PBS



## Basic Information

<b>Catalog Number:</b> 68253-8-PBS	<b>GenBank Accession Number:</b> BC030833	<b>Purification Method:</b> Protein G purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 4719	<b>CloneNo.:</b> 3C4E1
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P28331	
<b>Isotype:</b> IgG1	<b>Full Name:</b> NADH dehydrogenase (ubiquinone) Fe-S protein 1, 75kDa (NADH- coenzyme Q reductase)	
<b>Immunogen Catalog Number:</b> AG3135	<b>Calculated MW:</b> 727 aa, 79 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50476-4, NDUFS1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68253-8-PBS. Detection antibody: 68253-7-PBS. Standard: Ag3135. Range: 0.098-100 ng/mL.