

For Research Use Only

POLR3D Monoclonal antibody, PBS Only



Catalog Number: 68263-1-PBS

Featured Product

Basic Information

Catalog Number:

68263-1-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7646

GenBank Accession Number:

BC004484

GeneID (NCBI):

661

UNIPROT ID:

P05423

Full Name:

polymerase (RNA) III (DNA directed)

polypeptide D, 44kDa

Calculated MW:

44 kDa

Observed MW:

44 kDa

Purification Method:

Protein G purification

CloneNo.:

2F5G4

Applications

Tested Applications:

Indirect ELISA, WB

Species Specificity:

rabbit, human

Background Information

DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Polymerase (RNA) III (DNA directed) polypeptide D (POLR3D) plays a key role in sensing and limiting infection by intracellular bacteria and DNA viruses and acts as nuclear and cytosolic DNA sensor involved in innate immune response. POLR3D can sense non-self dsDNA that serves as template for transcription into dsRNA.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68263-1-Ig (POLR3D antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68263-1-PBS in a different storage buffer formulation.



WB result of POLR3D antibody (68263-1-Ig; 1:1000; incubated at room temperature for 1.5 hours) with sh-Control and sh-POLR3D transfected HEK-293 cells. This data was developed using the same antibody clone with 68263-1-PBS in a different storage buffer formulation.



1x10⁶ HeLa cells were intracellularly stained with 0.8 ug POLR3D Monoclonal antibody (68263-1-Ig; Clone:2F5G4) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1)(red), or 0.8 ug Mouse IgG1 isotype control Mouse McAb (66360-1-Ig; Clone: 1F8D3) (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011). This data was developed using the same antibody clone with 68263-1-PBS in a different storage buffer