

For Research Use Only

DHFR Monoclonal antibody, PBS Only



Catalog Number: 68274-1-PBS

Basic Information

Catalog Number: 68274-1-PBS	GenBank Accession Number: BC000192	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 1719	CloneNo.: 4C11E9
Source: Mouse	UNIPROT ID: P00374	
Isotype: IgG1	Full Name: dihydrofolate reductase	
Immunogen Catalog Number: AG7318	Calculated MW: 21 kDa	
	Observed MW: 21 kDa	

Applications

Tested Applications:
Indirect ELISA, WB

Species Specificity:
mouse, human

Background Information

DHFR (dihydrofolate reductase) catalyzes the subsequent NADPH-dependent reduction of the H₂folate produced by TS to regenerate the tetrahydrofolate (H₄folate) necessary for continued dTMP synthesis and transfer of 1-carbon units (PMID:3461439). It belongs to the dihydrofolate reductase family and can exist as a homodimer (PMID:2248959). DHFR is widely expressed in fetal and adult human brain and whole blood, expression is higher in adult brain than fetal brain (PMID:21310276).

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68274-1-Ig (DHFR antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 68274-1-PBS in a different storage buffer formulation.



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 68274-1-Ig (DHFR antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68274-1-PBS in a different storage buffer formulation.



1x10⁶ MCF-7 cells were intracellularly stained with 0.4 ug DHFR Monoclonal antibody (68274-1-Ig, Clone:4C11E9) and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1)(red), or 0.4 ug Mouse IgG1 isotype control Mouse McAb (66360-1-Ig, Clone: 1F8D3) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C). This data was developed using the same antibody clone with 68274-1-PBS in a

