

For Research Use Only

DHFR Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68274-7-PBS



Basic Information

Catalog Number: 68274-7-PBS	GenBank Accession Number: BC000192	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1719	CloneNo.: 4E6A5
Source: Mouse	UNIPROT ID: P00374	
Isotype: IgG1	Full Name: dihydrofolate reductase	
Immunogen Catalog Number: AG7318	Calculated MW: 21 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50686-3, DHFR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68274-6-PBS. Detection antibody: 68274-7-PBS. Standard:Ag7318. Range: 0.098-100 ng/mL.



Cytometric bead array standard curve of MP50686-4, DHFR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68274-8-PBS. Detection antibody: 68274-7-PBS. Standard:Ag7318. Range: 0.098-25 ng/mL.