

For Research Use Only

LEPREL2 Monoclonal antibody, PBS Only



Catalog Number: 68299-1-PBS

Basic Information

Catalog Number: 68299-1-PBS	GenBank Accession Number: BC017217	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 10536	CloneNo.: 1F6F3
Source: Mouse	UNIPROT ID: Q8IVL6	
Isotype: IgG1	Full Name: leprecan-like 2	
Immunogen Catalog Number: AG8850	Calculated MW: 736 aa, 82 kDa	
	Observed MW: 82 kDa	

Applications

Tested Applications:
WB, Indirect ELISA, IHC

Species Specificity:
Human, Rat, Mouse

Background Information

P3H3 (prolyl 3-hydroxylase 3, also known as LEPREL2) is a member of the leprecan family of proteins, which also include P3H1, P3H2, CRTAP and SC56. Collagen prolyl hydroxylases are required for proper collagen biosynthesis, folding, and assembly. P3H3 is thought to reside in the endoplasmic reticulum and has prolyl 3-hydroxylase activity catalyzing the post-translational formation of 3-hydroxyproline in -Xaa-Pro-Gly-sequences in collagens, especially types IV and V.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 68299-1-Ig (LEPREL2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68299-1-PBS in a different storage buffer formulation.



rat testis tissue were subjected to SDS PAGE followed by western blot with 68299-1-Ig (LEPREL2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68299-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded mouse ovary tissue slide using 68299-1-Ig (LEPREL2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 68299-1-PBS in a different storage buffer formulation.