

For Research Use Only

LEPREL2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68299-5-PBS



Basic Information

Catalog Number: 68299-5-PBS	GenBank Accession Number: BC017217	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 10536	CloneNo.: 1A1B4
Source: Mouse	UNIPROT ID: Q8IVL6	
Isotype: IgG1	Full Name: leprecan-like 2	
Immunogen Catalog Number: AG8850	Calculated MW: 736 aa, 82 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50849-2, LEPREL2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68299-4-PBS. Detection antibody: 68299-5-PBS. Standard: Ag8850. Range: 0.098-100 ng/mL.