

For Research Use Only

cTnT Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68300-3-PBS

Basic Information

Catalog Number: 68300-3-PBS	GenBank Accession Number: BC002653	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 7139	CloneNo.: 2C3E5
Source: Mouse	UNIPROT ID: P45379	
Isotype: IgG1	Full Name: troponin T type 2 (cardiac)	
Immunogen Catalog Number: AG7716	Calculated MW: 36 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50338-2, cTnT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68300-3-PBS. Detection antibody: 68300-2-PBS. Standard:Ag7716. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50338-3, cTnT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68300-3-PBS. Detection antibody: 68300-5-PBS. Standard:Ag7716. Range: 0.195-100 ng/mL.