

For Research Use Only

DLAT Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68303-5-PBS



Basic Information

Catalog Number:

68303-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4226

GenBank Accession Number:

BC039084

GeneID (NCBI):

1737

UNIPROT ID:

P10515

Full Name:

dihydrolipoamide S-acetyltransferase

Calculated MW:

647 aa, 70 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1B5D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50549-2, DLAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68303-4-PBS. Detection antibody: 68303-5-PBS. Standard:Ag4226. Range: 0.391-100 ng/mL.