

For Research Use Only

DLAT Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68303-8-PBS



Basic Information

| | | |
|-------------------------------------|--|--|
| Catalog Number: 68303-8-PBS | GenBank Accession Number: BC039084 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 1737 | CloneNo.: 1C12A5 |
| Source: Mouse | UNIPROT ID: P10515 | |
| Isotype: IgG1 | Full Name: dihydrolipoamide S-acetyltransferase | |
| Immunogen Catalog Number: AG4226 | Calculated MW: 647 aa, 70 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50549-4, DLAT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68303-8-PBS. Detection antibody: 68303-9-PBS. Standard: Ag4226. Range: 0.098-100 ng/mL.