

For Research Use Only

ATP5F1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68304-2-PBS



Basic Information

Catalog Number: 68304-2-PBS	GenBank Accession Number: BC005366	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 515	CloneNo.: 3D7E5
Source: Mouse	UNIPROT ID: P24539	
Isotype: IgG1	Full Name: ATP synthase, H ⁺ transporting, mitochondrial FO complex, subunit B1	
Immunogen Catalog Number: AG8571	Calculated MW: 256 aa, 29 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50484-1, ATP5F1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68304-2-PBS. Detection antibody: 68304-3-PBS. Standard:null. Range: 3.125-100 ng/mL.