

For Research Use Only

MED23 Monoclonal antibody, PBS Only



Catalog Number: 68307-1-PBS

Basic Information

Catalog Number: 68307-1-PBS	GenBank Accession Number: BC060759	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 9439	CloneNo.: 1H4G4
Source: Mouse	UNIPROT ID: Q9ULK4	
Isotype: IgG1	Full Name: mediator complex subunit 23	
Immunogen Catalog Number: AG28483	Calculated MW: 156 kDa	
	Observed MW: 150 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human, Mouse, Rat

Background Information

MED23 is a component of the Mediator complex, it is a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery. Mediator is recruited to promoters by direct interactions with regulatory proteins and serves as a scaffold for the assembly of a functional pre-initiation complex with RNA polymerase II and the general transcription factors.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 68307-1-Ig (MED23 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68307-1-PBS in a different storage buffer formulation.



4T1 cells were subjected to SDS PAGE followed by western blot with 68307-1-Ig (MED23 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 68307-1-PBS in a different storage buffer formulation.