

For Research Use Only

PNPT1 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 68309-4-PBS



Basic Information

Catalog Number:

68309-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6290

GenBank Accession Number:

BC053660

GeneID (NCBI):

87178

UNIPROT ID:

Q8TCS8

Full Name:

polyribonucleotide
nucleotidyltransferase 1

Calculated MW:

86 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3D6A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50567-2, PNPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68309-2-PBS. Detection antibody: 68309-4-PBS. Standard:Ag6290. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50567-3, PNPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68309-5-PBS. Detection antibody: 68309-4-PBS. Standard:Ag6290. Range: 0.781-100 ng/mL.