

For Research Use Only

# PNPT1 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68309-5-PBS



## Basic Information

Catalog Number:

68309-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG6290

GenBank Accession Number:

BC053660

GeneID (NCBI):

87178

UNIPROT ID:

Q8TCS8

Full Name:

polyribonucleotide  
nucleotidyltransferase 1

Calculated MW:

86 kDa

Purification Method:

Protein G purification

CloneNo.:

3F3D11

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50567-3, PNPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68309-5-PBS. Detection antibody: 68309-4-PBS. Standard:Ag6290. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP50567-4, PNPT1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68309-5-PBS. Detection antibody: 68309-3-PBS. Standard:Ag6290. Range: 0.098-6.25 ng/mL.