

For Research Use Only

# DDT Monoclonal antibody, PBS Only (Capture)

Catalog Number: 68314-5-PBS



## Basic Information

Catalog Number: 68314-5-PBS	GenBank Accession Number: BC015508	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1652	CloneNo.: 3B9E12
Source: Mouse	UNIPROT ID: P30046	
Isotype: IgG2a	Full Name: D-dopachrome tautomerase	
Immunogen Catalog Number: AG3065	Calculated MW: 118 aa, 13 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50531-2, DDT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68314-5-PBS. Detection antibody: 68314-4-PBS. Standard: Ag3065. Range: 0.781-100 ng/mL.