

For Research Use Only

NFIC Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68318-2-PBS

Basic Information

Catalog Number:

68318-2-PBS

Size:

1mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG9472

GenBank Accession Number:

BC012120

GeneID (NCBI):

4782

UNIPROT ID:

P08651

Full Name:

nuclear factor I/C (CCAAT-binding
transcription factor)

Calculated MW:

428aa, 48 kDa; 508aa, 57 kDa

Purification Method:

Protein A purification

CloneNo.:

1F6D10

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50142-1, NFIC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68318-2-PBS. Detection antibody: 68318-3-PBS. Standard: Ag9472. Range: 1.563-100 ng/mL.