

For Research Use Only

NFIC Monoclonal antibody, PBS Only (Detector)



Catalog Number: 68318-3-PBS

Basic Information

Catalog Number: 68318-3-PBS	GenBank Accession Number: BC012120	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4782	CloneNo.: 2F7G3
Source: Mouse	UNIPROT ID: P08651	
Isotype: IgG2a	Full Name: nuclear factor I/C (CCAAT-binding transcription factor)	
Immunogen Catalog Number: AG9472	Calculated MW: 428aa,48 kDa; 508aa,57 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50142-1, NFIC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68318-2-PBS. Detection antibody: 68318-3-PBS. Standard: Ag9472. Range: 1.563-100 ng/mL.