

For Research Use Only

MTHFD1L Monoclonal antibody, PBS Only (Capture)



Catalog Number: 68321-2-PBS

Basic Information

Catalog Number:

68321-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG9061

GenBank Accession Number:

BC017477

GeneID (NCBI):

25902

UNIPROT ID:

Q6UB35

Full Name:

methylenetetrahydrofolate
dehydrogenase (NADP+ dependent) 1-
like

Calculated MW:

978 aa, 106 kDa

Purification Method:

Protein G purification

CloneNo.:

2A7B1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50271-1, MTHFD1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68321-2-PBS. Detection antibody: 68321-3-PBS. Standard: Ag9061. Range: 0.195-100 ng/mL